



Prices and offers valid 7/1/18 through 9/30/18

SUMMERTIME SAVINGS

ROTO-THERM™ INCUBATED TUBE ROTATOR

- Digitally programmed to maintain chamber temperature (Amb +5 to 60°)
- Gentle fixed speed mixing at 24 RPM
- Small footprint 11.6" x 8.5" x 9"
- Includes common tube holders (0.5mL through 50mL)
- Optional tube holders are available



Cat. #	Description	LIST	SALE
H2020	Roto-Therm™ Incubated Rotator with tube holders	\$499.00	\$399.00

PROPETTE™ XL PIPETTE CONTROLLER

- Long-life rechargeable lithium battery for up to 10 hours continuous use
- Variable electronic speed adjustment
- Easy-fit silicone nosepiece accepts pipette sizes from 1mL to 100mL
- Low pressure fingertip control buttons
- Ergonomic body design, fits comfortably in either the right or left hand
- Comes with charger, benchtop stand & extra filter



Cat. #	Description	LIST	SALE
P6080	Propette Pipette Controller w/charger & US type plug	\$279.00	\$199.00

MC-24 TOUCH HIGH SPEED MICROCENTRIFUGE

- Full color touch screen
- mySpin™ one button program recall
- Combi Rotor accepts 24 tubes and two PCR strip tubes
- Variable speed to 13,500rpm / 16,800xg



Cat. #	Description	LIST	SALE
C2417	MC-24 Touch, Microcentrifuge 24x1.5ml / PCR Combi rotor	\$1,995.00	\$1,450.00

PCRPAL® PERSONAL THERMAL CYCLER

- Full featured stand-alone cycler, no PC or tablet needed
- Intuitive programming with color display and real time graphing
- Available with a 24 place 0.2mL or a 20 place 0.5mL block
- Economically priced, ideal for school and research labs
- Small footprint; Only 7" x 9" x 13"D



Cat. #	Description	LIST	SALE
ATS4005	PCRPal Personal Thermal with 25 x 0.2ml block	\$995.00	\$875.60

PROPETTE™ LE ADJUSTABLE PIPETTES

- Designed for accuracy and comfort
- Features an adjustable volume with Roto-Lock™
- Superior accuracy and outstanding user comfort
- Full volume range from 0.1µl to 10mL



Cat. #	Description	LIST	SALE
P5200-2U	ProPette LE single channel adjustable Pipet, 0.2-2µl, Ea	\$219.00	\$175.00
P5200-10U	ProPette LE single channel adjustable Pipet, 1-10µl, Ea	\$219.00	\$175.00
P5200-20U	ProPette LE single channel adjustable Pipet, 2-20µl, Ea	\$219.00	\$175.00
P5200-50U	ProPette LE single channel adjustable Pipet, 5-50µl, Ea	\$219.00	\$175.00
P5200-100U	ProPette LE single channel adjustable Pipet, 10-100µl, Ea	\$219.00	\$175.00
P5200-200U	ProPette LE single channel adjustable Pipet, 20-200µl, Ea	\$219.00	\$175.00
P5200-300U	ProPette LE single channel adjustable Pipet, 30-300µl, Ea	\$219.00	\$175.00
P5200-1M	ProPette LE single channel adjustable Pipet, 100-1000µl, Ea	\$219.00	\$175.00

100mm PETRI DISHES

- Extremely flat, high clarity polystyrene
- Stackable for easy storage and handling
- Sterilized, available in single, 2, 3 or 4 compartments



Cat. #	Description	LIST	SALE
704005	100 x 15mm Petri Dish, Sterile, Stackable, 20/bag, 500/cs	\$ 81.14	\$60.00
M090201	Petri Dish, 150x15mm, PS, STRL, Slippable, 100/cs	\$ 88.00	\$66.00
M090501	Deep Dish, 100x20mm, STRL, Slippable, 20/sleeve, 400/cs	\$122.00	\$90.50

CELL CULTURE DISHES

- High clarity polystyrene
- Flat transparent surface for distortion-free observation
- Vacuum plasma TC treatment; excellent cell adherence
- PA/PE packaging; Stackable
- Vents for gas exchange; Sterile; Non-pyrogenic



Cat. #	Description	LIST	SALE
704001	Cell Culture Dish, 100x20mm, TC, Sterile, 20/bag, 300/CS	\$ 98.53	\$ 73.89
715001	Cell Culture Dish, 150 x 20mm, TC, Sterile, 5/pk, 100/CS	\$139.10	\$104.32
705001	Cell Culture Dish, 60 x 15mm, TC, Sterile, 20/pk, 500/CS	\$107.32	\$ 80.49
706001	Cell Culture Dish, 35 x 10mm, TC, Sterile, 20/pk, 500/CS	\$ 92.74	\$ 69.55

SURESTRAIN™ PREMIUM CELL STRAINERS

- Advanced design for faster flow through
- Superior ventilation prevents clogging and spills
- Sterile, RNase, DNase, pyrogen and cytotoxin free
- Packed in individual blister packs



Cat. #	Description	LIST	SALE
M854040	SureStrain™ Premium Cell Strainers, 40µm, Blue, 50/pk	\$ 99.00	\$74.25
M854100	SureStrain™ Premium Cell Strainers, 100µm, Yellow, 50/pk	\$ 99.00	\$74.25
M854070	SureStrain™ Premium Cell Strainers, 70µm, White, 50/pk	\$ 99.00	\$74.25

READY-STRAIN™ CELL STRAINING KIT

- Includes 50mL centrifuge tube with attached strainer, cap
- Available in 40, 70 and 100µm pore sizes
- Individually wrapped, sterile, 50 per pack

Cat. #	Description	LIST	SALE
M855040	Ready Strain Cell Straining Kit, 40µm, Str1, 50mL tube, 50/pk	\$139.00	\$104.25
M855070	Ready Strain Cell Straining Kit, 70µm, Str1, 50mL tube, 50/pk	\$139.00	\$104.25
M855100	Ready Strain Cell Straining Kit, 100µm, Str1, 50mL tube, 50/pk	\$139.00	\$104.25



Laboratory Products Sales Inc.

A NY State Certified Woman-Owned Small Business

1665 Buffalo Road • Rochester, N.Y. 14624 • 585-247-4720

TOLL FREE: 1-800-388-0166 • FAX: 585-247-6686 • www.LPSinc.com