



LAB ESSENTIALS SALE



CellPro™
Premium Labware



Vacuum Filtration Systems and Bottle Top Filters

- Includes the filter funnel and the polystyrene flasks
- Individually wrapped, 12 per case and gamma sterilized
- Durable PES membrane for low protein binding and fast filtering
- Available with 0.22 or 0.45µm filter membranes



Cat#	Description	List	SALE PRICE
V05022	CellPro™ Vacuum Filtration System, 50mL, 0.22µm, case/12	\$99.95	\$50.75
V05045	CellPro™ Vacuum Filtration System, 50mL, 0.45µm, case/12	\$99.95	\$50.75
V15022	CellPro™ Vacuum Filtration System, 150mL, 0.22µm, case/12	\$99.95	\$50.75
V15045	CellPro™ Vacuum Filtration System, 150mL, 0.45µm, case/12	\$99.95	\$50.75
V25022	CellPro™ Vacuum Filtration System, 250mL, 0.22µm, case/12	\$129.95	\$62.35
V25045	CellPro™ Vacuum Filtration System, 250mL, 0.45µm, case/12	\$129.95	\$62.35
V50022	CellPro™ Vacuum Filtration System, 500mL, 0.22µm, case/12	\$149.95	\$89.00
V50045	CellPro™ Vacuum Filtration System, 500mL, 0.45µm, case/12	\$149.95	\$89.00
V100022	CellPro™ Vacuum Filtration System, 1L, 0.22µm, case/12	\$199.95	\$145.00
V100045	CellPro™ Vacuum Filtration System, 1L, 0.45µm, case/12	\$199.95	\$145.00
VH05022	CellPro™ Vacuum Filtration BTF, 50mL, 0.22µm, case/24	\$99.95	\$79.96
VH05045	CellPro™ Vacuum Filtration BTF, 50mL, 0.45µm, case/24	\$99.95	\$79.96
VH15022	CellPro™ Vacuum Filtration BTF, 150mL, 0.22µm, case/24	\$99.95	\$79.96
VH15045	CellPro™ Vacuum Filtration BTF, 150mL, 0.45µm, case/24	\$99.95	\$79.96
VH25022	CellPro™ Vacuum Filtration BTF, 250mL, 0.22µm, case/24	\$129.95	\$103.96
VH25045	CellPro™ Vacuum Filtration BTF, 250mL, 0.45µm, case/24	\$129.95	\$103.96
VH50022	CellPro™ Vacuum Filtration BTF, 500mL, 0.22µm, case/24	\$149.95	\$119.96
VH50045	CellPro™ Vacuum Filtration BTF, 500mL, 0.45µm, case/24	\$149.95	\$119.96
VH100022	CellPro™ Vacuum Filtration BTF, 1L, 0.22µm, case/24	\$199.95	\$159.96
VH100045	CellPro™ Vacuum Filtration BTF, 1L, 0.45µm, case/24	\$199.95	\$159.96

PurePoint™ Precision Pipet Tips

- Low Binding Fine Tip Design Insures the Most Accurate Pipetting
- Works Well With Most Pipet Brands Including Gilson® and Eppendorf®
- Extremely Straight Tip For Use With 8 or 12 Channel Pipettors
- Certified RNase, DNase, Human DNA Free
- Available in 96 Tip Sterile Hinged Boxes or Deck Reloads



Cat. #	Description	Unit	List	SALE
L1RT101	PurePoint Deck Reloads 10µL (xl)	10x96	\$23.00	\$17.00
L1RT120	PurePoint Deck Reloads 200µL	10x96	\$23.00	\$17.00
L1RT130	PurePoint Deck Reloads 300µL	10x96	\$23.00	\$17.00
L1RT125	PurePoint Deck Reloads 1250µL	8x96	\$23.00	\$17.00
Starter Pack (Sterile Tips in Hinged Boxes)				
L1ST101	PurePoint Deck Reloads 10µL (xl)	10x96	\$35.95	\$21.00
L1ST120	PurePoint Deck Reloads 200µL	10x96	\$35.95	\$21.00
L1ST130	PurePoint Deck Reloads 300µL	10x96	\$35.95	\$21.00
L1ST125	PurePoint Deck Reloads 1250µL	8x96	\$35.95	\$21.00

SURESEAL S™ STERILE MICROCENTRIFUGE TUBES

- Available in 1.5mL, 2.0mL, and 0.5mL volumes
- Sterile packed in resealable self-standing bags
- Graduated BioClear™ tubes for maximum clarity
- High Speed (28,000xg) and Autoclavable (121°C)
- Includes 10 each Stop-Pops™ Cap Locks in each bag



Cat. #	Description	LIST	SALE PRICE
L220001	SureSeal S, MCT, 1.5mL, Grad STRL, Clear, 500/bg	\$12.90	\$9.00
L220021	SureSeal S, MCT, 2.0mL, Grad STRL, Clear, 500/bg	\$14.00	\$10.00
L220051	SureSeal S, MCT, 1.5mL, Grad STRL, Clear, 50/pk, 500 /CS	\$17.00	\$13.00
L220071	SureSeal S, MCT, 0.5mL, Grad STRL, Clear, 500/bg	\$12.90	\$9.00



Laboratory Products Sales Inc.

A NY State Certified Woman-Owned Small Business

1665 Buffalo Road • Rochester, N.Y. 14624 • 585-247-4720

TOLL FREE: 1-800-388-0166 • FAX: 585-247-6686 • www.LPSinc.com